

# Intensity Gain from VULCAN Guide 9

April 2013



A U.S. Department of Energy Multilaboratory Project

SPALLATION NEUTRON SOURCE  
SNS  
Argonne National Laboratory • Brookhaven National Laboratory • Thomas Jefferson National Accelerator Facility • Lawrence Berkeley National Laboratory • Los Alamos National Laboratory • Oak Ridge National Laboratory

This report was prepared as an account of work sponsored by an agency of the United States government. Neither the United States government nor any agency thereof, nor any of their employees, makes any warranty, express or implied, or assumes any legal liability or responsibility for the accuracy, completeness, or usefulness of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights. Reference herein to any specific commercial product, process, or service by trade name, trademark, manufacturer, or otherwise, does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States government or any agency thereof. The views and opinions of authors expressed herein do not necessarily state or reflect those of the United States government or any agency thereof.

## **Intensity Gain from VULCAN Guide 9**

E. B. Iverson

A. D. Stoica

W. Lu

Date Published: 26 April 2013

Prepared by Oak Ridge National Laboratory  
for  
UT-BATTELLE, LLC  
under contract DE-AC05-00OR22725  
for the  
U.S. DEPARTMENT OF ENERGY

Intensity Gain from VULCAN Guide 9		Page 1 of 3
106100200-TR0???-R00	Author: E. B. Iverson	Reviewer: F. X. Gallmeier

## 1. INTRODUCTION

VULCAN is the SNS engineering diffractometer, [1] and has been operating successfully within the user program for some time. Guide segment 9 is being commissioned at present; this report documents the observed increase in flux attainable at the sample position with this final guide segment in place.

## 2. INTENSITY MEASUREMENTS

We used a calibrated beam monitor [2] and the default instrument chopper settings to characterize the spectral current passing through a small pinhole (2.32 mm in diameter) located at the nominal sample position. We performed these measurements before and after the installation of guide segment 9, with the upstream guide configuration in both the High-Intensity and High-Resolution modes. The measured intensities appear in Figure 1, and the data are embedded in the electronic version of this document.

The ratio before and after Guide 9 installation is shown in Figure 2, and again the data are attached to the electronic version of this document. Guide 9 has minimal effect in the High-Resolution mode (as intended). In High-Intensity mode, guide 9 increases the measured flux by around 25% at 1 Å and nearly 70% at 5 Å, consistent with design estimates.<sup>1</sup>

### 2.1 HISTORICAL COMPARISON

Although this study was carried out to document the performance gain due to guide segment 9, we also obtain absolute intensity numbers we can compare to earlier measurements. The initial commissioning measurements [1] quote a peak spectral intensity (measured in the High-Intensity configuration) of  $2.7 \times 10^8 \text{ n Å}^{-1} \text{ MW}^{-1} \text{ s}^{-1}$ . This value integrates over the entirety of a beam spot of area  $20.6 \text{ cm}^2$ . It was not reported in [1] but was measured at the same time (using image plates and narrow-band chopper settings) that the peak-to-average ratio of the spatial distribution across that beam spot was 3.5, giving a peak flux of  $4.6 \times 10^7 \text{ n cm}^{-2} \text{ Å}^{-1} \text{ MW}^{-1} \text{ s}^{-1}$ . Our new measurement, as shown in Figure 1(a) above, gives a peak flux of  $4.1 \times 10^7 \text{ n cm}^{-2} \text{ Å}^{-1} \text{ MW}^{-1} \text{ s}^{-1}$ . While somewhat smaller now, this reduction is consistent with other measurements showing that the intensity per unit accelerator power has some dependence on accelerator power—that is, intensity is slightly less than linear with accelerator power.

## 3. PEAK SHAPE MEASUREMENTS

[Put in measured increase in peak shape width?]

## 4. SUMMARY


Guide segment 9 on the VULCAN instrument at SNS has been installed and characterized. It results in a flux gain of nearly 70% at long wavelengths and more than 20% at 1 Å in the VULCAN High-Intensity configuration. The absolute intensity measured in the VULCAN sample position now is comparable to the value measured at instrument commissioning, and peaks in wavelength at  $4.1 \times 10^7 \text{ n cm}^{-2} \text{ Å}^{-1} \text{ MW}^{-1} \text{ s}^{-1}$  when measured at 0.8 MW across a pinhole of 2.3 mm diameter.

## REFERENCES


- [1] X.-L. Wang, T. M. Holden, A. D. Stoica, K. An, H. D. Skorpenske, A. B. Jones, G. Q. Rennich, and E. B. Iverson, “First results from the VULCAN diffractometer at the SNS,” *Materials Science Forum*, vol. 652, pp. 105–110, 2010.
- [2] E. B. Iverson, B. J. Micklich, D. V. Baxter, R. G. Cooper, P. D. Ferguson, D. W. Freeman, F. X. Gallmeier, S. E. Hammons, C. M. Lavelle, and I. Popova, “Neutronic measurements to commission

---

<sup>1</sup>Need reference.



(a) High-Intensity configuration.



(b) High-Resolution configuration.

Figure 1: Spectral fluxes measured with and without guide segment 9.




Figure 2: Intensity gain measured for guide segment 9.

the SNS,” in *Proceedings of ICANS XVII, the Seventeenth Meeting of the International Collaboration on Advanced Neutron Sources*, 2005.